



JETi[®]

Thrombectomy System

Powerful
THROMBECTOMY for the
Toughest Peripheral Clot





Powerful Thrombectomy for the Toughest Peripheral Clot

- ✓ Rapid, aggressive clot removal
- ✓ Gentle to the vessel wall
- ✓ Designed for a single session procedure with or without adjunctive thrombolytics



Saline jet inside catheter macerates clot without clinically significant hemolysis

Clot Detect™ – advanced clot sensing provides flow feedback to reduce blood loss

HyperPulse™ – optional thrombolytic infusion mode¹



Federal law restricts the JETi device to sale by or on the order of a licensed physician

¹ HyperPulse™ is CE marked, not cleared by the FDA
MM0006-OUS, Rev. A



Arterial JETi Case




Thrombosed femoral stent




Post JETi thrombectomy, vessel fully patent


Venous JETi Case



Thrombosed Iliac vein



Post JETi thrombectomy and balloon maceration



Post stent, vein fully patent

JETi® Thrombectomy System

Description	REF Number	Guide Wire Size
Sterile Disposables		
JETi 8Fr Catheter (100 cm)	WV8100	0.014" - 0.035"
JETi Pump Set	WV06PS01	
Pressure Monitoring	K11-01392	
Suction Tubing	28-0206	
Non-Sterile Disposables		
Canister 850 cc	HCS7850	
Bacteria Filter	01-90-3928	
Capital Equipment		
JETi Saline Drive Unit	WV0008	
Vacuum pump	01-23-0405	



USA Toll Free: +1-844-796-8324

International: +1-949-422-1222

Email: support@walkvascular.com

www.JETi8.com